

**LUDLUM MODEL 2360
SCALER/RATEMETER DATA LOGGER**

March 2020

**Serial Number 342290 and Succeeding
Serial Numbers**

**LUDLUM MODEL 2360
SCALER/RATEMETER DATA LOGGER**

March 2020

**Serial Number 342290 and Succeeding
Serial Numbers**



LUDLUM MEASUREMENTS, INC
501 OAK STREET, P.O. BOX 810
SWEETWATER, TEXAS 79556
325-235-5494, FAX: 325-235-4672

STATEMENT OF WARRANTY

Ludlum Measurements, Inc. warrants the products covered in this manual to be free of defects due to workmanship, material, and design for a period of twelve months from the date of delivery. The calibration of a product is warranted to be within its specified accuracy limits at the time of shipment. In the event of instrument failure, notify Ludlum Measurements to determine if repair, recalibration, or replacement is required.

This warranty excludes the replacement of photomultiplier tubes, G-M and proportional tubes, and scintillation crystals which are broken due to excessive physical abuse or used for purposes other than intended.


There are no warranties, express or implied, including without limitation any implied warranty of merchantability or fitness, which extend beyond the description of the face there of. If the product does not perform as warranted herein, purchaser's sole remedy shall be repair or replacement, at the option of Ludlum Measurements. In no event will Ludlum Measurements be liable for damages, lost revenue, lost wages, or any other incidental or consequential damages, arising from the purchase, use, or inability to use product.

RETURN OF GOODS TO MANUFACTURER

If equipment needs to be returned to Ludlum Measurements, Inc. for repair or calibration, please send to the address below. All shipments should include documentation containing return shipping address, customer name, telephone number, description of service requested, and all other necessary information. Your cooperation will expedite the return of your equipment.

**LUDLUM MEASUREMENTS, INC.
ATTN: REPAIR DEPARTMENT
501 OAK STREET
SWEETWATER, TX 79556**

| REV # | ALTERATIONS | DATE | BY |
|-------|-----------------|---------|-----|
| 1 | VALID | 8/1/96 | RSS |
| 2 | UPDATED ARTWORK | 5/15/06 | CMC |




| | | |
|--|----------------------|---------------------|
| DMN DATE CMC 5-15-06 | CHK DATE | APP DATE 5-15-06 |
| PART NUM: 4390-159 | SCALE: FULL OR OTHER | |
| TITLE M 2360 ALPHA/BETA DATA LOGGER | | |
| LUDLUM MEASUREMENTS, INC. 501 ONE STREET SWEETWATER, TEXAS 75086 | SERIES 390 | SHEET 157 |

Table of Contents

| | |
|---|-----------------|
| <i>Introduction</i> | <i>1</i> |
| <i>Getting Started</i> | <i>2</i> |
| Unpacking and Repacking | 2-1 |
| Battery Installation | 2-1 |
| Connecting a Detector to the Instrument | 2-2 |
| Internal Switches | 2-2 |
| Battery Test | 2-2 |
| Operating the Instrument | 2-3 |
| Principle of Operation | 2-3 |
| <i>Specifications</i> | <i>3</i> |
| <i>Description of Controls and Functions</i> | <i>4</i> |
| Operator Controls | 4-1 |
| Internal Controls | 4-3 |
| <i>Safety Considerations</i> | <i>5</i> |
| Environmental Conditions for Normal Use | 5-1 |
| Cleaning Instructions and Precautions | 5-1 |
| Warning Markings and Symbols | 5-1 |
| <i>Calibration and Maintenance</i> | <i>6</i> |
| Calibration | 6-1 |
| Establishing an Operating Point | 6-1 |
| Meter Calibration | 6-2 |
| Detector Overload Calibration | 6-3 |
| Maintenance | 6-4 |
| Recalibration | 6-4 |
| Batteries | 6-5 |
| <i>RS-232 Interface</i> | <i>7</i> |
| Communicating With the Model 2360 | 7-1 |
| Commands | 7-1 |
| Read Commands | 7-1 |

| | |
|--|------------------|
| Set Commands | 7-2 |
| Command Descriptions | 7-2 |
| Model 2360 Interface Software | 7-6 |
| Functions | 7-7 |
| Main Screen | 7-8 |
| Auto Dump Data (Setup) | 7-9 |
| Model 2360 Interface Software Sample Printouts | 7-10 |
| Installation of the Model 2360 Interface Software | 7-12 |
| Removal of the Model 2360 Interface Software | 7-13 |
| Connecting to a Computer | 7-13 |
| Software License Agreement | 7-14 |
| | |
| <i>Technical Theory of Operation</i> | <i>8</i> |
| Amplifier/Power Supply Board | 8-1 |
| Processor Board | 8-3 |
| | |
| <i>Troubleshooting</i> | <i>9</i> |
| Troubleshooting Electronics that Utilize Proportional & Scintillator Type Detectors | 9-1 |
| | |
| <i>Recycling</i> | <i>10</i> |
| | |
| <i>Parts List</i> | <i>11</i> |
| Model 2360 Scaler/Ratemeter | 11-1 |
| Amplifier/Power Supply Board, Drawing 390 × 400 | 11-1 |
| Processor Board, Drawing 390 × 395 | 11-4 |
| LCD Display Board, Drawing 390 × 372 | 11-6 |
| Calibration Board, Drawing 390 × 176 | 11-6 |
| BCD Board, Drawing 261 × 107 | 11-6 |
| Wiring Diagram, Drawing 390 × 179 | 11-7 |
| | |
| <i>Drawings and Diagrams</i> | <i>12</i> |



Introduction

The Ludlum Model 2360 Scaler/Ratemeter Data Logger is an easy-to-use survey instrument incorporating the best of the analog and digital worlds. Able to measure both alpha and beta radiation levels simultaneously, the Model 2360 presents the data in scaler (digital display) mode or ratemeter (analog meter) mode. The Model 2360 also logs up to 550 data points, consisting of sample number, date/time stamp, sample measurements (both alpha and beta, cpm if ratemeter measurement, or counts per count time if logging scaler count), scaler count time, “S” or “R,” identifying whether scaler count or ratemeter was logged, and location identifier (10 characters).

Each alpha-beta data point is logged by simply pressing the button located in the instrument handle. The appropriate scaler or ratemeter measurement is taken, and both alpha and beta readings are logged into non-volatile memory. An internal dipswitch allows for logging of the scaler reading, ratemeter reading, both scaler and ratemeter readings, and NO logging. The location identifier can be input by a PC or other RS-232 device prior to logging a sample.

There are also 6 lines (15 characters each) of header information that can be stored at the beginning of the non-volatile memory. The header lines can contain such information as the user name, survey name, serial numbers, etc. The “CALIBRATION DUE DATE” can also be stored in non-volatile memory, which will disable the instrument when the internal clock date reaches the stored date.

The Model 2360 has a long arc length (6.0 cm {2.4 in.}) meter face that normally reads from 0 to 500 cpm. The main rotary switch allows for multiplication ranges of $\times 1$, $\times 10$, $\times 100$, and $\times 1000$. The Model 2360 has a three-position toggle switch on the front panel to switch between displaying alpha, beta, or alpha + beta levels.

Alpha and beta pulses have different audio tones, so that the user can discriminate between the two. A two-position momentary action switch allows either a meter reset, or displays the detector HV onto the meter. The following six alarms may be set via the RS-232 port:

Alpha Ratemeter (0-999999 cpm)

Beta Ratemeter (0-999999 cpm)

Alpha + Beta Ratemeter (0-999999 cpm)

Alpha Scaler (0-999999)

Beta Scaler (0-999999)


Alpha + Beta Scaler (0-999999)

These alarm points, when exceeded, light the meter face LED, marked ALARM, and activate the audio speaker. The ratemeter alarms are non-latching (will cease when the radiation level drops below the alarm point). The scaler alarm will continue until the RESET is pressed, or the next scaler count is started.

The digital display is a full 6-digit liquid crystal display (LCD) display, which is direct-driven for good viewing and a wide temperature response from -4 to 60 °C (40 to 140 °F). The digital display displays the scaler count, and prior to each sample logging, displays the current sample number. The display also has an arrow symbol for counting overflow and two colons that indicate that a scaler count is in progress.

The Model 2360 communicates through an RS-232 port located on the instrument can. The computer interface software is supplied on CD and includes the following functions:

- downloading of header and logged data into an ASCII file
- setting and reading of instrument parameters/header information
- setting of the "CALIBRATION DUE DATE"
- setting of the internal real-time clock/calendar
- setting of the user-definable scaler time
- setting of location code and alarm points
- clearing of logged memory

A dark gray square with the word "Section" in white at the top and a large white number "2" in the center.

Getting Started

Unpacking and Repacking

Remove the calibration certificate and place it in a secure location. Remove the instrument and accessories (batteries, cable, etc.) and ensure that all of the items listed on the packing list are in the carton. Check individual item serial numbers and ensure calibration certificates match. The Model 2360 serial number is located on the front panel below the battery compartment. Most Ludlum Measurements, Inc. detectors have a label on the base or body of the detector for model and serial number identification.

Important!


If multiple shipments are received, ensure that the detectors and instruments are not interchanged. Each instrument is calibrated to specific detectors, and is therefore not interchangeable.

To return an instrument for repair or calibration, provide sufficient packing material to prevent damage during shipment. Also provide appropriate warning labels to ensure careful handling

Every returned instrument must be accompanied by an **Instrument Return Form**, which can be downloaded from the Ludlum website at www.ludlums.com. Find the form by clicking the “Support” tab and selecting “Repair and Calibration” from the drop-down menu. Then choose the appropriate Repair and Calibration division where you will find a link to the form.

Battery Installation

Ensure the Model 2360 range selector switch is in the “OFF” position. Open the battery lid by pushing down and turning the quarter-turn thumbscrew



counterclockwise a quarter of a turn. Install two "D" size batteries in the compartment.

Note the (+) and (-) marks inside the battery door. Match the battery polarity to these marks. Close the battery box lid, push down, and turn the quarter-turn thumb screw clockwise a quarter of a turn.


Note:

The center post of a flashlight battery is positive. The batteries are placed in the battery compartment in opposite directions.

Connecting a Detector to the Instrument

Caution!

The detector operating voltage (HV) is supplied to the detector via the detector input connector. A mild electric shock may occur if you make contact with the center pin of the input connector. Switch the Model 2360 to the "OFF" position before connecting or disconnecting the cable or detector.



Connect one end of a detector cable to the detector by firmly pushing the connectors together while twisting clockwise $\frac{1}{4}$ turn. Repeat the process in the same manner with the other end of the cable and the instrument.

Internal Switches

Release the can latches and remove the cover from the 2360, taking care not to damage the speaker wires. Using a ball-point pen, set the switches for the desired AUDIO division and TONE described on Page 4-3. Replace the cover and fasten the latches.

Battery Test

The batteries should be checked each time the instrument is turned on. Move the range selector switch to the BAT position. Ensure that the meter needle deflects to the battery check portion on the meter scale. If the meter does not respond, check to see if the batteries have been correctly installed. Replace the batteries if necessary.

Turn the range selector switch from the off position to the $\times 1000$ position and verify that the meter needle is driven full-scale for approximately two seconds and then returns to zero. Also ensure both alarm LEDs on the front panel turn on briefly. The LCD should go through an initialization sequence displaying "88:8.8:8.8," then the current sample number, and finally "0."

Operating the Instrument

Connect a detector to the instrument if you have not already done so. Obtain a meter reading from a check source or calibrated source, if available. Verify that the reading falls within the expected range. Remove the source.

If a radiation source is available, increase the meter count to exceed the alarm threshold. Both the appropriate alarm lamp and audio alarm signal should activate.

Depress the RESET toggle switch. The meter needle should drive to zero and the alarm circuit should de-energize, shutting off both visual and audible alarms.

Proceed with use.


Principle of Operation

The Model 2360 is to be used in combination with an alpha/beta scintillation or proportional detector. The pulse height differential between the alpha and beta pulses from the scintillation or proportional detector is discriminated by the Model 2360.

The detected alpha count is displayed by selecting the " α " position on the three-position " $\alpha/\alpha+\beta/\beta$ " toggle switch. Likewise, the sum of the alpha and beta counts or the beta counts only are displayed by selecting the appropriate " $\alpha + \beta$ " or " β " position. Multiply the cpm reading on the analog ratemeter by the range multiplier position. When using the LCD and preset count time interval, the counts are accumulated in each of the three channels during the count cycle. The count cycle is started by depressing the push-button located in the end of the carrying handle. If a scaler or ratemeter alarm is activated, the ALARM LED will light and the audio will sound a steady tone.

The RESET toggle switch position resets the meter pointer to zero and deactivates any current alarm. The detector operating voltage is displayed on

the meter dial by selecting the READ HV position. The high-voltage scale ranges from 0 to 2 kV. The OL (overload) lamp, located in the lower, right corner of the meter dial, indicates that the detector is saturated either by a puncture in the detector face or an exposure to a radiation field above the counting capability of the instrument. The analog meter deflects full scale when the OL lamp is illuminated.

A dark gray square with the word "Section" in white at the top and a large white number "3" in the center.

Specifications

Compatible Detectors: proportional and dual phosphor scintillation detectors; common Models: 43-1-1, 43-2-2, 43-20, 43-68, 43-89 and 43-93

Data Logger: capable of logging up to 550 individual data points into non-volatile memory with the following identifiers for each point:

- alpha and beta sample counts
- sample number
- date/time stamp
- scaler count time
- 10-character location identifier

Range: four linear range multiples of $\times 1$, $\times 10$, $\times 100$, and $\times 1000$; used in combination with the 0-500 CPM meter dial providing an overall range of 0-500 kcpm

Thresholds: internal control allows for adjustment from -2 to -15 mV for beta and -40 to -700 mV for alpha

Window: internal control allows for adjustment from the beta threshold up to the alpha threshold setting (beta only)

Audio: built-in unimorph speaker with volume control (greater than 60 dB at 0.61 m {2 ft}, full volume); headset jack located on the instrument can

Audio Divide: selectable dual or individual click-per-event for alpha and beta counts with divisions of 1, 10, 100, or 1000 events per click (beta only)

Alarm points: six separate alarm points, set through the RS-232 port, activating the alarm audio tone and lighting the ALARM LED. The six alarm points can be set from 0 to 99999. Ratemeter alarms are non-latching, while scaler alarms are not. (A latched alarm requires the RESET button be pressed in order for the alarm to clear) The six alarms are:

- Alpha Ratemeter
- Beta Ratemeter
- Alpha + Beta Ratemeter
- Alpha Scaler
- Beta Scaler
- Alpha + Beta Scaler

High Voltage: recessed front-panel potentiometer; adjustable from 200 to 2000 Vdc

Linearity: within 10% of true value for the analog CPM meter; within 2% for the LCD

Response Time: $\times 1 = 7$ seconds, $\times 10 = 2$ seconds, $\times 100 = 2$ seconds, $\times 1000 = 2$ seconds; all response times measured from 10-90% of full scale

Analog Meter: rugged 1 milliamp (mA), with pivot-and-jewel suspension and 6.0 cm (2.4 in.) arc length

Connector: series "C" standard, others available

Temperature Range: -20 to 50 °C (-4 to 122 °F); may be certified for operation from -40 to 65 °C (-40 to 150 °F)

Maximum Relative Humidity: less than 95% (non-condensing)


Power: two standard "D" size batteries housed in an externally accessible sealed compartment

Battery Life: typically 250 hours of operation with a fresh set of alkaline "D" cell batteries

Size: 16.5 x 8.9 x 21.6 cm (6.5 x 3.5 x 8.5 in.) (H x W x L)

Weight: 1.6 kg (3.5 lb), including batteries

Finish: drawn-and-cast aluminum, with beige powder coating

A dark gray square containing the word "Section" in a small, white, sans-serif font at the top, and a large, white, bold number "4" in the center.

Description of Controls and Functions

Operator Controls

OFF/BAT/X1000/X100/X10/X1 Switch (or Range Selector Switch): a six-position rotary switch to select the analog meter range multipliers and check the battery status. When switched to the "BAT" position, the meter pointer should deflect above the left vertical mark on the "BAT OK" line. Moving the range selector switch to one of the range multiplier positions ($\times 01$, $\times 10$, $\times 100$, $\times 1000$) provides the operator with an overall range of 0-500 kcpm. Multiply the scale reading by the multiplier to determine the actual reading.

During initial instrument turn-on, the meter will be driven full scale for about two seconds and then return to zero. The LCD will show "88:8.8:8.8," display the current sample number, and then "0." If the count button is pressed while the display is showing all 8s then the following numbers are shown: firmware number, date, time, PC defined scaler time, alpha ratemeter alarm, alpha + beta rate meter alarm, alpha scaler alarm, alpha + beta scaler alarm, current sample number, and then "0."

Liquid Crystal Display (LCD): six-digit display that shows the scaler count for the selected channel or the sum of both. It also indicates when a count is in progress by displaying two colons; the colons are turned off when the count is complete. If the counter exceeds 999999, an arrow in the upper-left corner of the display turns on in order to indicate an overflow; the counter then rolls over to zero and continues counting.

VOL: Turning this control clockwise will increase the speaker volume, and counterclockwise will decrease the volume.

Note:

The volume should be turned down, when not required, to reduce battery drain.

$\alpha/\alpha+\beta/\beta$ Switch: A three-position toggle switch used to select the sum of alpha and beta count channels ($\alpha+B$), alpha count only (α), or beta count only (B), for display. This switch affects both the analog meter and digital display. The rate meter channels are active regardless of the switch position and will continue to function when the channel is not selected for display. This ability allows the operator to view each channel separately or together by simply selecting the appropriate switch position.

Headphone Jack: 0.32 cm (0.13 in. or one-eighth inch) phone jack on instrument can for the connection of external headphones. Inserting a phone plug into the jack disconnects the external unimorph speaker from the audio circuitry. Use 0.32 cm (0.13 in. or one-eighth inch) mating plug-LMI part # 21-9653.

RS-232 Port: located on the instrument can and allows for connection of the instrument to a PC for dumping data and setting up parameters

RESET/HV: a dual-position momentary toggle switch which provides readout of the detector high voltage when the HV position is selected. When the RESET position is selected, this switch provides a rapid means of driving the analog ratemeter to zero and resetting the meter after an alarm condition. Use the 0-2 kV meter scale for high-voltage readings.

MIN 0.1, 0.5, 1, 2, 5, 10, 60, PC Count Time Select Switch: an eight-position rotary switch used to select scaler count times in minutes. When set to PC, the scaler will use the user-defined count time. The count cycle is initiated by depressing the push-button switch in the carrying handle.

Log Push-button (located in the carrying handle): When depressed, the current sample number is displayed, the counter is reset to zero, and the timer is started. If selected, the instrument will also log the current reading. The colons on the display turn on and stay on until the count has expired.

Remove the front panel calibration cover to access the following calibration potentiometer:

HV: a multi-turn potentiometer used to adjust the detector high voltage from 200 to 2000 Vdc

Internal Controls

Remove the instrument cover (can) to access the following dipswitches on SW1.

RECYCLE: A one-pole DIP switch (#1) used to select recycle scaler mode. When placed in the ON position, the instrument will start a count cycle. If the instrument is set to log the sample, it will be saved and a new count will start again. The ratio is selected from the following table:

| SWITCH 1 | RECYCLE MODE |
|---------------------------|-------------------------------|
| ON | ON |
| OFF | OFF |

TONE: a one-pole DIP switch (#2) used to select tone discrimination between alpha and beta count channels. When in the DUAL mode, alpha and beta pulse tones will be audible in all selector switch positions (i.e. if in the β -only position and beta radiation is detected, the beta tones will be heard in addition to the alpha tones).

When the SNGL tone position is selected, both alpha and beta pulse tones can be heard in the “ $\alpha+\beta$ ” selection, but alpha pulses are the only audible tones in the α -channel selection, and beta pulse tones are the only audible tones in the β -channel selection.

| SWITCH 2 | TONE MODE |
|---------------------------|----------------------------|
| ON | DUAL |
| OFF | SINGLE |

AUDIO Divide: A two-pole DIP switch (#3 & #4) used to select audio divisions of 1, 10, 100, and the beta audio subtract mode.

Note:

The AUDIO divide function only affects the lower-frequency beta tones. The higher-frequency alpha click-per-event will be unaffected by the divide-by selection.

The ratio is selected from the following table:

| SWITCH | | DIVIDE BY |
|--------|-----|-----------------------------|
| 3 | 4 | RATIO |
| ON | ON | 1 |
| ON | OFF | 10 |
| OFF | ON | 100 |
| OFF | OFF | Beta audio Subtract Mode |

LOG SAMPLE: A two-pole DIP switch (#5 & #6) used to select logging mode. When both switches are OFF, the instrument will not log samples but will output the ratemeter readings to the RS-232 port every two seconds in ASCII format. When set to log both scaler and ratemeter, two samples are actually saved. The logging mode is selected from the following table.

| SWITCH | | LOGGING |
|--------|-----|--|
| 5 | 6 | MODE |
| ON | ON | Log ratemeter and scaler |
| ON | OFF | Log ratemeter |
| OFF | ON | Log scaler |
| OFF | OFF | No sample saved. Output ratemeter every 2 seconds. |

The following controls are utilized during calibration only and should only be altered by a qualified calibration technician.

MTR: a multi-turn potentiometer used to calibrate the meter to the CPM reading

AT: a multi-turn potentiometer used to vary the alpha pulse discriminator from 40 to 700 millivolts (mV)

BW: a multi-turn potentiometer used to vary the upper beta pulse discriminator from BT setting to AT setting


BT: a multi-turn potentiometer used to vary the lower beta pulse discriminator from 2 to 15 mV

OL: a multi-turn potentiometer that provides a means to vary the detector current overload set point

LIM: a multi-turn potentiometer used to set the maximum HV limit to 2000 Vdc

HV: a multi-turn potentiometer used to adjust the high-voltage test reading (0 to 2 kV scale) to correspond with the actual high-voltage output. The HV switch must be depressed during adjustment.

LB: a multi-turn potentiometer used to adjust the minimum battery voltage level (2.2 Vdc), corresponding to the low-battery indication on the meter dial. The BAT switch position must be selected during adjustment.

A dark gray square with the word "Section" in white at the top and a large white number "5" in the center.

Safety Considerations

Environmental Conditions for Normal Use

Indoor or outdoor use

No maximum altitude

Temperature range of -20 to 50 °C (-4 to 122 °F); may be certified for operation from -40 to 65°C (-40 to 150 °F)

Maximum relative humidity of less than 95% (non-condensing)

Pollution Degree 3 (as defined by IEC 664) (Occurs when conductive pollution or dry nonconductive pollution becomes conductive due to condensation. This is typical of industrial or construction sites.)

Cleaning Instructions and Precautions

The Model 2360 Scaler/Ratemeter may be cleaned externally with a damp cloth, using only water as the wetting agent. Do not immerse the instrument in any liquid. Observe the following precautions when cleaning:

1. Turn the instrument range selector switch to the OFF position.
2. Allow the instrument to sit for one minute before cleaning.

Warning Markings and Symbols

Caution!

The operator or responsible body is cautioned that the protection provided by the equipment may be impaired if the equipment is used in a manner not specified by Ludlum Measurements, Inc.

Caution!

Verify instrument voltage input rating before connecting to a power converter. If the wrong power converter is used, the instrument and/or power converter could be damaged.

The Model 2360 Scaler/Ratemeter is marked with the following symbols:

CAUTION (per ISO 3864, No. B.3.1) – designates hazardous live voltage and risk of electric shock. During normal use, internal components are hazardous live. This instrument must be isolated or disconnected from the hazardous live voltage before accessing the internal components. This symbol appears on the front panel. **Note the following precautions:**

Warning!

The operator is strongly cautioned to take the following precautions to avoid contact with internal hazardous live parts that are accessible using a tool:

1. Turn the instrument power OFF and disconnect the power cord.
2. Allow the instrument to sit for one minute before accessing internal components.




CAUTION, RISK OF ELECTRIC SHOCK (per ISO 3864, No. B.3.6) – designates a terminal (connector) that allows connection to a voltage exceeding 1 kV. Contact with the subject connector while the instrument is on or shortly after turning off may result in electric shock. This symbol appears on the front panel.



The “**crossed-out wheellie bin**” symbol notifies the consumer that the product is not to be mixed with unsorted municipal waste when discarding; each material must be separated. The symbol is placed on the battery lid. See section 10, “Recycling,” for further information.



The “CE” mark is used to identify this instrument as being acceptable for use within the European Union. This is located on the front panel.

A dark gray square containing the word "Section" in a small, white, sans-serif font at the top, and a large, white, bold number "6" in the center.

Calibration and Maintenance

Calibration

Note:

Local procedures may supersede the following.

ESTABLISHING AN OPERATING POINT

The detector operating parameters are established by adjusting the detector operating voltage (HV), alpha threshold, and beta window to find an optimum efficiency for the alpha/beta scintillator or proportional detector.

The threshold and window parameters can be adjusted to optimize alpha/beta count discrimination, count efficiency, and minimize "cross talk" between channels. Refer to the specific detector operation manual or calibration certificate for the suggested threshold and window settings. Once the thresholds and window settings are established, an operating voltage-versus-count rate plot should be performed for both alpha and beta count channels with alpha and beta particle emission sources.

The following procedure is an example of determining the operating voltage for an alpha/beta scintillation or proportional detector:

- Connect a Ludlum Model 500 Pulser or equivalent to the Model 2360.
- Switch the 2360 to the β position. Adjust the beta threshold (BT) for 3.5 mV and the window (BW) for 30 mV. The pulser counts should be detected on the 2360 ratemeter above 3.5 ± 1 mV and should shut off above 30 mV.
- Switch the channel selector switch to the α position. Adjust the pulser for 120 mV pulse output and vary the AT control until counts are detected on the ratemeter.

- Depress the HV switch and adjust the HV potentiometer for 0.4 to 0.5 kV on the 0-2.0 kV scale. Connect the scintillator and switch to the β only position. Place an alpha source on the detector face.
- Slowly increase the HV potentiometer to observe an increase, then decrease, and increase again in count as the HV is increased. Decrease the HV until the ratemeter is in the "dip" of the observed count rate-versus-HV plot just performed. Depress the HV switch and note the HV setting.
- Plot a HV versus count rate plateau in 25 V increments, 50 V each side of the HV reading found in the above step (i.e., HV setting for count "dip" in the above step = 675 V, start the plot at 625 V and increase in 25 V steps until 725 V is reached). Plot alpha source, beta source, and background counts for both α and β channel positions.
- Find the optimum operating voltage from the plot, which gives the greatest alpha and beta source efficiency while maintaining no greater than the maximum acceptable level of cross talk between channels.
- Select the desired count channel display, and proceed to use instrument.

METER CALIBRATION

A Ludlum Model 500 Pulser or equivalent is required. If the pulser does not have high-voltage readout, use a high-impedance voltmeter with at least 1000 megohm input resistance to measure the detector voltage.

Ensure that the meter movement has proper mechanical zero. The adjustment is on the front of the meter bezel. It must be adjusted to "0" with the ON/OFF selector switch in the OFF position.

Connect the Model 500 Pulser to the Model 2360 with the appropriate cable. Rotate the Model 2360 range multiplier switch to the $\times 100$ position. Select the $\alpha + \beta$ channel position.

Adjust the Pulser for 40,000 cpm and adjust the pulse amplitude to at least twice the beta threshold level (i.e.; beta threshold = 3.5 mV, adjust the pulser to 7-10 mV).

Remove the instrument cover and adjust the MTR potentiometer until the meter reads 400 cpm. Adjust the pulser to 10,000 cpm and ensure ratemeter reads $100 \pm 10\%$. Decade the pulser and Model 2360 range multiplier switch

to check meter linearity on the $\times 1000$, $\times 10$, and $\times 1$ positions. Linearity should be within 10% of each reading.

Set the LCD count time for one minute. Adjust the pulser count rate to 40 kcpm. Depress the count button, and when the count cycle is complete, confirm that LCD reads within 2% of the incoming count rate.

Adjust the BT, BW, and AT controls for the appropriate set points as described in the following subsection.

Connect a high-impedance, high-voltage meter (may use the Model 500 Pulser if equipped with an HV meter) to the detector input connector and adjust the HV control for a reading of 1000 Vdc on the voltmeter.

Depress the RES/HV switch to the HV position and adjust the HV potentiometer located on the circuit board for a reading of 1.0 kV on the meter dial. Adjust the HV output from 500 to 1500 Vdc and confirm that the 2360 HV meter corresponds to the external voltmeter within 10% of each reading.

Rotate the range multiplier switch to the OFF position. Remove the batteries from the battery compartment and connect a DC power supply to the two screw terminals located at the rear of the battery compartment. The positive power supply lead should connect to the terminal with the red wire, and the negative lead to the terminal with the black wire.

Adjust the power supply for 2.2 Vdc and switch the 2360 to the BAT position. Adjust the LB potentiometer to align the meter needle with the low-battery mark on the meter dial (vertical line to the left of BAT OK).

Replace the instrument can and proceed with use.

DETECTOR OVERLOAD CALIBRATION

Note:

The detector operating voltage (HV) must be determined and set before the OL (overload) adjustment is performed. If the detector operating voltage is re-adjusted, the overload adjustment must be re-adjusted.

Adjust the OL control to the maximum counterclockwise position.

Note:

Detector saturation is when the meter response no longer increases with increasing radiation field intensity.

For alpha/beta scintillators, expose the detector photomultiplier tube (PMT) to a small light leak by loosening the detector window. Some scintillation detectors incorporate a screw in the detector body, which when removed, will simulate a detector face puncture or light leak. The meter should start to decrease toward zero as light saturates the scintillation material.

Expose just enough light to where the meter starts to decrease. Adjust the OL control until the overload LED just begins to flicker on the meter dial. The ratemeter should deflect above full-meter scale at this point.

Re-seal the detector window and expose the detector to a radiation source that will drive the meter near full scale. Confirm that the LED does not turn on and the meter remains on scale.

Maintenance

Instrument maintenance consists of keeping the instrument clean and periodically checking the batteries and the calibration. The Model 2360 instrument may be cleaned with a damp cloth (using only water as the wetting agent). Do not immerse instrument in any liquid. Observe the following precautions when cleaning:

1. Turn the instrument OFF and remove the batteries.
2. Allow the instrument to sit for one minute before accessing internal components.

RECALIBRATION

Recalibration should be accomplished after maintenance or adjustments have been performed on the instrument. Recalibration is not normally required following instrument cleaning, battery replacement, or cable replacement.

Note:

Ludlum Measurements, Inc. recommends recalibration at intervals no greater than one year. Check the appropriate regulations to determine required recalibration intervals.


Ludlum Measurements offers a full-service repair and calibration department. We not only repair and calibrate our own instruments, but most other manufacturers' instruments as well. Calibration procedures are available upon request for customers who choose to calibrate their own instruments.

BATTERIES

The batteries should be removed any time the instrument is placed into storage. Battery leakage may cause corrosion on the battery contacts, which must be scraped off and/or washed using a paste solution made from baking soda and water. Use a spanner wrench to unscrew the battery contact insulators, exposing the internal contacts and battery springs. Removal of the handle will facilitate access to these contacts.

Note:

Never store the instrument over 30 days without removing the batteries. Although this instrument will operate at very high ambient temperatures, battery seal failure may occur at temperatures as low as 37 °C (100 °F).

A dark gray square with the word "Section" in white at the top and a large white number "7" in the center.

RS-232 Interface

Communicating with the Model 2360

The Model 2360 has an RS-232 serial port that can be connected to a PC or other RS-232 device. The proper communication settings are (2400,8,N,1):

2400 Baud
8 data bits
No parity
1 stop bit

Commands

The following table shows all the commands used to communicate with the Model 2360. All commands must be in uppercase letters. The data can be in either uppercase or lowercase. All set commands should return "OK" + a carriage return and line feed.

READ COMMANDS

| | |
|-----|-------------------------------------|
| RA | Read sample number |
| RC | Read calendar date |
| RD | Read date |
| RHx | Read header |
| RL | Read current location |
| RP | Read user (PC) time |
| RR | Read α / β ratemeter |
| RS | Read samples |
| R1 | Read α ratemeter alarm |
| R2 | Read β ratemeter alarm |
| R3 | Read $\alpha+\beta$ ratemeter alarm |
| R4 | Read α scaler alarm |
| R5 | Read β scaler alarm |
| R6 | Read $\alpha+\beta$ scaler alarm |

SET COMMANDS

| | |
|--------------------|------------------------------------|
| SCmmddyyyy | Set calibration date |
| SDmmddyyyy | Set current date |
| SHxxxxxxxxxxxxxxxx | Set Header |
| SLxxxxxxxxxxx | Set current location |
| SMx | Sets ratemeter dumping |
| SPxxx.x | Set user (PC) time |
| SR | Send reset samples |
| SThhmm | Set current time |
| S1xxxxxx | Set α ratemeter alarm |
| S2xxxxxx | Set β ratemeter alarm |
| S3xxxxxx | Set $\alpha+\beta$ ratemeter alarm |
| S4xxxxxx | Set α scaler alarm |
| S5xxxxxx | Set β scaler alarm |
| S6xxxxxx | Set $\alpha+\beta$ scaler alarm |

COMMAND DESCRIPTIONS**RA**

This command reads the current sample number. The output is six characters, including a [CR] and [LF]. The format is:

```
0001 [CR] [LF]
```

RC

This command reads the calibration date. During power-up, the Model 2360 checks the current date against this date. If the current date is past the calibration date, the message "OUTCAL" is displayed and the Model 2360 is disabled until the calibration date is changed ahead of the current date. The output is 10 characters including a [CR] and [LF]. The format is:

```
MMDDYYYY  
12251996 [CR] [LF]
```

RD

This command reads the current date and time. The output is 21 characters including a [CR] and [LF]. The format is :

```
HH:MM:SS MM/DD/YYYY  
16:16:29 11/20/1996 [CR] [LF]
```

RHx

This command reads the specified header where “x” equals a number 1-6. The output is 17 characters, including a carriage return [CR] & line-feed [LF].

```
[15 characters]  
John Smith [5 SPACES] [CR] [LF]
```

RL

This command reads the current location. The output is 12 characters, including a [CR] and [LF]. The format is:

```
TABLE0001 [SPACE] [CR] [LF]
```

RP

This command reads the user-defined “PC” count time. This is the scaler count time when the Model 2360 count time switch is on the “PC” position. The output is 7 characters, including a [CR] and [LF]. The format is:

```
012.5 [CR] [LF]
```

RR

This command reads the current alpha and beta ratemeter reading. When dip-switch 5 and 6 are set to **OFF**, the Model 2360 sends this message every 2 seconds. The alpha Ratemeter reading is first, followed by the beta ratemeter reading. The output is 15 characters, including a [CR] and [LF]. The format is:

```
Alpha   Beta  
000003 002305 [CR] [LF]
```

RS

This command will return all logged samples from memory. A “\$” signifies the end of samples. The maximum number of samples stored is 550. The format is:

```
0001 11/18/96 14:50:05 000020 000450 R 001.0 CHKSRC  
0002 11/18/96 14:50:07 000015 000390 S 001.0 TABLE0001  
0003 11/18/96 16:49:49 000040 001400 R 000.1 TABLE0002  
$
```

R1, R2, R3, R4, R5, R6

These commands read the alarm set points.

R1 = Alpha Ratemeter.
 R2 = Beta Ratemeter.
 R3 = Alpha + Beta Ratemeter.
 R4 = Alpha Scaler.
 R5 = Beta Scaler.
 R6 = Alpha + Beta Scaler.

The output is eight characters, including a [CR] and [LF]. The format is:

```
000500 [CR] [LF]
```

SCmmdyyy

This command sets the calibration due date. The date is entered in Month-Day-Year (MMDDYYYY) format. During power-up, the Model 2360 checks to see if the current date is past the calibration due date. If it is, then the Model 2360 displays “OUTCAL” and is disabled until this command is issued to set the calibration date ahead of the current date. The length of the command is 12 characters, including a [CR] and [LF]. The format is:

```
SC11201997 [CR] [LF]
```

Sdmmddyyyz

This command sets the current date. The date is entered in Month-Day-Year format. The PCF8593 clock/calendar chip uses a counter from 0 to 3 to represent the year. The variable “z” must correspond to the following table. The length of the command is 13 characters, including a [CR] and [LF].

| Year ending in | Year Code |
|----------------|-----------|
| 04 | 0 |
| 05 | 1 |
| 06 | 2 |
| 07 | 3 |
| 08 | 0 |
| 09 | 1 |
| 10 | 2 |
| 11 | 3 |

For example, the command to set the date to March 20, 2006 is:

```
SD032020062 [CR] [LF]
```

SHxyyyyyyyyyyyyyyy

This command sets the specified header. The variable “x” can be any number between 1 and 6. The variable “y” must be 15 characters. If the value is less than 15 characters, it must be padded with spaces. The length of the command is 20 characters, including a [CR] and [LF]. The format is:

```
SH1JOHN SMITH[5 SPACES] [CR] [LF]
```

SLxxxxxxxxxx

This command sets the current location that will be saved with any subsequent logged samples. The location can be up to 10 characters in length and must be padded with spaces if less than 10 characters. The length of the command is 14 characters, including a [CR] and [LF]. The format is:

```
SLTABLE0007[SPACE] [CR] [LF]
```

SMx

This command disables or enables the automatic dumping of the ratemeter when the Model 2360 is set not to log samples (dip-switch 5 and 6 both OFF). Specifying SM0 will disable the ratemeter dumping until the unit is turned off or the command SM1 is issued.

SPxxx.x

This command sets the user-defined count time that is selectable by setting the count time switch on the Model 2360 to “PC”. The count time can be set from 000.1 minutes (6 seconds) to 546.1 minutes (32766 seconds). The length of the command is nine characters, including a [CR] and [LF]. The format is:

```
SP001.0 [CR] [LF]
```

SR

This command resets the sample number to one and clears all samples stored in memory. Use this command with caution. The length of the command is four characters, including a [CR] and [LF]. The format is:

WARNING!

This command will erase all logged samples from memory.

```
SR[CR] [LF]
```

SThhmm

This command sets the current time in 24-hour format. Twenty-four-hour time is as follows:

| | | | |
|----------|------|----------|------|
| 12:00 AM | 0000 | 12:00 PM | 1200 |
| 01:00 AM | 0100 | 01:00 PM | 1300 |
| 02:00 AM | 0200 | 02:00 PM | 1400 |
| 03:00 AM | 0300 | 03:00 PM | 1500 |
| 04:00 AM | 0400 | 04:00 PM | 1600 |
| 05:00 AM | 0500 | 05:00 PM | 1700 |
| 06:00 AM | 0600 | 06:00 PM | 1800 |
| 07:00 AM | 0700 | 07:00 PM | 1900 |
| 08:00 AM | 0800 | 08:00 PM | 2000 |
| 09:00 AM | 0900 | 09:00 PM | 2100 |
| 10:00 AM | 1000 | 10:00 PM | 2200 |
| 11:00 AM | 1100 | 11:00 PM | 2300 |

The length of the command is eight, characters including a [CR] and [LF]. For example, the command to set the time to 1:00 pm is:

ST1300[CR][LF]

S1, S2, S3, S4, S5, S6

These commands set the alarms for the alpha, beta, and alpha + beta ratemeter and also the alpha, beta, and alpha + beta scaler. If the alarm is set to 0, then the alarm is disabled. S1 = Alpha Ratemeter.

S2 = Beta Ratemeter.

S3 = Alpha + Beta Ratemeter.

S4 = Alpha Scaler.

S5 = Beta Scaler.

S6 = Alpha + Beta Scaler.

The length of the command is 10 characters, including a [CR] and [LF]. For example, to set the alpha ratemeter alarm to 500 the command is:

S1000500 [CR] [LF]

Model 2360 Interface Software

The Model 2360 Interface (LMI Part #:1370-039) is Windows-based and has a user-friendly interface, which allows the user to communicate with the Model 2360. The Model 2360 interface features automatically loading default values, and a Auto Dump Mode Display. When the program is

started, the user is prompted to either allow the software to find the Model 2360 or to manually specify a serial port. Once connected, the software will download all data from the Model 2360, which includes headers and logged data. The user is able to change any information and update the Model 2360, print hard copies, or save data to an ASCII file for later import into word processors, spreadsheets, or other applications.

FUNCTIONS

The Model 2360 Interface software has three main functions:

1. Read or update the parameters stored in the Model 2360.
2. Read, save, or delete the logged data stored in the Model 2360.
3. Collect and save real-time data at specified intervals of time.

Model 2360 Interface version 2.1

User Defined Settings

Header 1: John Q. Public [Set]
 Header 2: SN: 220859 [Set]
 Header 3: SN: PR200747 [Set]
 Header 4: Site: Bldg 1 [Set]
 Header 5: RM 008, S. Wall [Set]
 Header 6: Comment [Set]
 Current Location: Table 007 [Set]

Date Settings

2360 Date: 03/17/2006 [Set]
 Current [Set]
 2360 Time: 16:29 [Set]
 Calibration Due Date: 03/17/2007 [Set]

Ratemeter Settings

Alpha Alarm: 999999 [Set]
 Beta Alarm: 999999 [Set]
 A+B Alarm: 999999 [Set]

Scaler Settings

Alpha Alarm: 999999 [Set]
 Beta Alarm: 999999 [Set]
 A+B Alarm: 999999 [Set]

Scaler Count Time (minutes)

User (PC) Time: 0.1 [Set]

Log Data - All readings in Counts Per Minute (CPM)

Next Sample Number: 33

| Sample # | Date | Time | Alpha | Beta | S/R | Count Time | Location |
|----------|------------|-------------|--------|------|-----|------------|-----------|
| 1 | 03/17/2006 | 04:13:25 PM | 311480 | 0 | R | 0.1 | Table 007 |
| 2 | 03/17/2006 | 04:13:27 PM | 31147 | 0 | S | 0.1 | Table 007 |
| 3 | 03/17/2006 | 04:13:39 PM | 311440 | 0 | R | 0.1 | Table 007 |
| 4 | 03/17/2006 | 04:13:41 PM | 31146 | 0 | S | 0.1 | Table 007 |
| 5 | 03/17/2006 | 04:13:53 PM | 311450 | 0 | R | 0.1 | Table 007 |
| 6 | 03/17/2006 | 04:13:55 PM | 31147 | 0 | S | 0.1 | Table 007 |
| 7 | 03/17/2006 | 04:14:07 PM | 311450 | 0 | R | 0.1 | Table 007 |
| 8 | 03/17/2006 | 04:14:09 PM | 31147 | 0 | S | 0.1 | Table 007 |

3/17/2006 4:52 PM

MAIN SCREEN

[**Set**]**]**—Clicking the **Set** button will save the parameter to the Model 2360. Each parameter has its own **Set** button.

[**Auto Dump**]**]**—Displays the Auto Dump Data screen, which allows for real time logging of data.

[**Clear Log Data**]**]**—Erases the logged sample data in the Model 2360.

Caution:

Data will be lost if this button is pressed and the data has not previously been saved.

[**Reload All Data**]**]**—Downloads all information from the Model 2360.


[**Load Defaults**]**]**—Reloads all header data to the original state as shipped from Ludlum Measurements, Inc.

[**Save Log Data**]**]**—Displays the "save file" prompt to allow the user to specify which drive and directory to save the Log File.

[**Save Log Data**]**]**—Displays the "save file" prompt to allow the user to specify which drive and directory to save the Log File. This option also saves the six header fields to the file.

[**Print**]**]**—Prints all parameters as well as logged data to the default printer. If there is no logged data, only the parameters will be printed.


[**Comm Setup**]**]**—Displays the Select Comm Port screen. Select "Automatic" to allow the software to scan all available serial ports to find the Model 2360. Select "Manual" and choose a specific serial port.



AUTO DUMP DATA (SETUP)

Note:

For auto dumping to work properly, the Model 2360 internal dip switches 5 and 6 must be set to the OFF position.



| Alpha | Beta | Date | Time |
|--------|------|------------|-------------|
| 311588 | 0 | 03/17/2006 | 04:59:07 PM |
| 311584 | 0 | 03/17/2006 | 04:59:09 PM |
| 311554 | 0 | 03/17/2006 | 04:59:11 PM |
| 311589 | 0 | 03/17/2006 | 04:59:13 PM |
| 311554 | 0 | 03/17/2006 | 04:59:15 PM |
| 311537 | 0 | 03/17/2006 | 04:59:17 PM |
| 311596 | 0 | 03/17/2006 | 04:59:19 PM |
| 311584 | 0 | 03/17/2006 | 04:59:21 PM |
| 311554 | 0 | 03/17/2006 | 04:59:23 PM |
| 311537 | 0 | 03/17/2006 | 04:59:25 PM |
| 311588 | 0 | 03/17/2006 | 04:59:27 PM |
| 311584 | 0 | 03/17/2006 | 04:59:29 PM |
| 311554 | 0 | 03/17/2006 | 04:59:31 PM |
| 311529 | 0 | 03/17/2006 | 04:59:33 PM |
| 311588 | 0 | 03/17/2006 | 04:59:35 PM |
| 311539 | 0 | 03/17/2006 | 04:59:37 PM |
| 311536 | 0 | 03/17/2006 | 04:59:39 PM |
| 311536 | 0 | 03/17/2006 | 04:59:41 PM |
| 311558 | 0 | 03/17/2006 | 04:59:43 PM |
| 311552 | 0 | 03/17/2006 | 04:59:45 PM |

All readings are in counts per minute.

Setup

2 Second Counts Second Count

60 Second Count

Count Down:

Readings

Averaged Discrete

Start Clear

Auto Scroll Grid

The Model 2360 internal dip-switch positions 5 and 6 must be switched to the OFF position to work properly

Save Log Data

Save Log Data with Headers

Print

Close

[Two-Second Counts]—The Model 2360 outputs the ratemeter reading every two seconds. This option will capture each two-second reading.

[User Defined Count]—This option will log the data based on a user-defined count time, which is adjustable from 2 to 9998 seconds. The count time specified here must be evenly divided by 2. The count time may be changed while a count is in progress and will take effect when the current count time is finished. When a user-defined count is in progress, the “Count Down” box will display the remaining count time.

[Readings]--When counting with a user-defined count time, the saved readings can be averaged or discrete. When the "Averaged" option is selected, the readings accumulated every two seconds during the count time are averaged. When the "Discrete" option is selected, the last two-second reading at the end of the count time is used. This option is only available when the user-defined count time option is selected.

[Start] [Stop]—Toggle from Start/Stop to start or stop the dumping of counts.

[Clear]—Clears data from the display box. This only clears the grid and does not affect the samples stored in the memory of the Model 2360.

[Auto Scroll Grid]—When checked, the grid automatically scrolls to keep the newest record visible.

[Save Log Data]—Save the data to a user-specified file and location. The data is saved in standard ASCII and is readable in almost any program.

[Save Log Data with Headers]—This button performs the same function as above, but also saves the six header fields.

[Print]—Sends data, including the displayed readings, to the default printer. If there is no data displayed, only the header information will be printed.

[Close]—Returns to the main menu area.

MODEL 2360 INTERFACE SOFTWARE SAMPLE PRINTOUTS (following pages)

Table 1: Logged data downloaded from the Model 2360 and saved to an ASCII file

Table 2: Auto dump data saved to an ASCII file

Table 3: Printout of logged data

Table 1:

Header 1: John Q. Public
 Header 2: SN: 220859
 Header 3: SN: PR200747
 Header 4: Site: Bldg 1
 Header 5: RM 008, S. Wall
 Header 6: Comment

S=Scaler, R=Rateometer

| Sample # | Date | Time | Alpha | Beta | S/R | Count Time | Location |
|----------|------------|-------------|--------|------|-----|------------|-----------|
| 1 | 03/17/2006 | 11:48:37 AM | 311092 | 0 | R | 0.1 | Table 007 |
| 2 | 03/17/2006 | 11:48:39 AM | 31112 | 0 | S | 0.1 | Table 007 |
| 3 | 03/17/2006 | 11:48:51 AM | 311144 | 0 | R | 0.1 | Table 007 |
| 4 | 03/17/2006 | 11:48:53 AM | 31113 | 0 | S | 0.1 | Table 007 |
| 5 | 03/17/2006 | 11:49:05 AM | 311127 | 0 | R | 0.1 | Table 007 |
| 6 | 03/17/2006 | 11:49:07 AM | 31113 | 0 | S | 0.1 | Table 007 |
| 7 | 03/17/2006 | 11:49:19 AM | 311143 | 0 | R | 0.1 | Table 007 |
| 8 | 03/17/2006 | 11:49:21 AM | 31113 | 0 | S | 0.1 | Table 007 |
| 9 | 03/17/2006 | 11:49:33 AM | 311129 | 0 | R | 0.1 | Table 007 |
| 10 | 03/17/2006 | 11:49:35 AM | 31113 | 0 | S | 0.1 | Table 007 |

Table 2:

Header 1: John Q. Public
 Header 2: SN: 220859
 Header 3: SN: PR200747
 Header 4: Site: Bldg 1
 Header 5: RM 008, S. Wall
 Header 6: Comment

| Alpha | Beta | Date | Time |
|--------|------|------------|-------------|
| 311348 | 0 | 03/17/2006 | 01:49:44 PM |
| 311344 | 0 | 03/17/2006 | 01:49:46 PM |
| 311299 | 0 | 03/17/2006 | 01:49:48 PM |
| 311348 | 0 | 03/17/2006 | 01:49:50 PM |
| 311314 | 0 | 03/17/2006 | 01:49:52 PM |
| 311312 | 0 | 03/17/2006 | 01:49:54 PM |
| 311297 | 0 | 03/17/2006 | 01:49:56 PM |
| 311296 | 0 | 03/17/2006 | 01:49:58 PM |
| 311348 | 0 | 03/17/2006 | 01:50:00 PM |
| 311314 | 0 | 03/17/2006 | 01:50:02 PM |

Table 3:

Model 2360 Log Data Date: 03/17/2006 Time: 01:44:49 PM Page 1

 Header 1: John Q. Public
 Header 2: SN: 220859
 Header 3: SN: PR200747
 Header 4: Site: Bldg 1
 Header 5: RM 008, S. Wall
 Header 6: Comment

Calibration Due Date: 03/15/2007
 Model 2360 Date: 03/17/2006
 Model 2360 Time: 01:37:21 PM

Logged Samples: 10

User PC Scaler Count Time: 0.1 minutes

Alpha Ratemeter Alarm Setpoint: 999999
 Beta Ratemeter Alarm Setpoint: 999999
 Alpha + Beta Ratemeter Alarm Setpoint: 999999

Alpha Scaler Alarm Setpoint: 999999
 Beta Scaler Alarm Setpoint: 999999
 Alpha + Beta Scaler Alarm Setpoint: 999999

S=Scaler, R=Ratemeter

| Sample # | Date | Time | Alpha | Beta | S/R | Count Time | Location |
|----------|------------|-------------|--------|------|-----|------------|-----------|
| 1 | 03/17/2006 | 11:48:37 AM | 311092 | 0 | R | 0.1 | Table 007 |
| 2 | 03/17/2006 | 11:48:39 AM | 31112 | 0 | S | 0.1 | Table 007 |
| 3 | 03/17/2006 | 11:48:51 AM | 311144 | 0 | R | 0.1 | Table 007 |
| 4 | 03/17/2006 | 11:48:53 AM | 31113 | 0 | S | 0.1 | Table 007 |
| 5 | 03/17/2006 | 11:49:05 AM | 311127 | 0 | R | 0.1 | Table 007 |
| 6 | 03/17/2006 | 11:49:07 AM | 31113 | 0 | S | 0.1 | Table 007 |
| 7 | 03/17/2006 | 11:49:19 AM | 311143 | 0 | R | 0.1 | Table 007 |
| 8 | 03/17/2006 | 11:49:21 AM | 31113 | 0 | S | 0.1 | Table 007 |
| 9 | 03/17/2006 | 11:49:33 AM | 311129 | 0 | R | 0.1 | Table 007 |
| 10 | 03/17/2006 | 11:49:35 AM | 31113 | 0 | S | 0.1 | Table 007 |

INSTALLATION OF THE MODEL 2360 INTERFACE SOFTWARE

User must comply with the software license agreement located at the end of this section (pages 7-14 thru 7-16). By installing this software you are consenting to be bound by this agreement. If you do not agree to all the terms of this agreement, do not install the product!

Insert the Model 2360 Interface software CD into the computer. The installation routine should start automatically. If it does not, click on the Start button, select "Run" and type in the following: "d:\setup.exe".

Replace the drive letter with the correct drive letter of the CD-ROM drive. Follow the onscreen prompts to install the software. When complete, the software should be installed in C:\Program Files\2360Win. A shortcut is created in the Start Menu under “Ludlum Measurements, Inc.”

REMOVAL OF THE MODEL 2360 INTERFACE SOFTWARE

To remove the program, start from the Windows Desktop. Click on Start Button, Settings, Control Panel, and then Add/Remove Programs. Find LMI 2360 Interface from the list and highlight, then click on the add/remove button to start uninstall.

Setup will prompt to ensure removal of program LMI 2360 Interface is acceptable.

Setup will also present a prompt asking if you wish to keep, remove, remove none, or remove all. Files can be removed since they are installed in C:\Program Files\Model 2360 and should not cause any problems.

Uninstall screen appears and the program removal will be complete.

CONNECTING TO A COMPUTER

Using the supplied cable, connect the end with the female connector to the Model 2360 and the other end to a free COM port on your computer. The pin-outs of the cables are as follows:

9-pin cable

| <u>Model 2360</u> | <u>Computer</u> |
|-------------------|-----------------|
| 2 | 2 |
| 3 | 3 |
| 5 | 5 |
| 7 | 7 |
| 8 | 8 |

25-pin cable

| <u>Model 2360</u> | <u>Computer</u> |
|-------------------|-----------------|
| 2 | 3 |
| 3 | 2 |
| 5 | 7 |
| 7 | 4 |
| 8 | 5 |



LUDLUM MEASUREMENTS, INC.

501 OAK ST., P.O. BOX 810

SWEETWATER, TX 79556

325/235-5494 FAX: 325/235-4672

Software License Agreement

Rev. (number) 1.0

Written by (or Revised by):

Date: 20 Jan 06

Approved by:

Richard Troala
Ed Em

Date: 20 Jan 06

BY INSTALLING THIS SOFTWARE, YOU ARE CONSENTING TO BE BOUND BY THIS AGREEMENT. IF YOU DO NOT AGREE TO ALL OF THE TERMS OF THIS AGREEMENT, DO NOT INSTALL THE PRODUCT.

Single User License Grant: Ludlum Measurements, Inc. ("Ludlum") and its suppliers grant to Customer ("Customer") a nonexclusive and nontransferable license to use the Ludlum software ("Software") in object code form solely on a single central processing unit owned or leased by Customer or otherwise embedded in equipment provided by Ludlum.

Customer may make one (1) archival copy of the Software provided Customer affixes to such copy all copyright, confidentiality, and proprietary notices that appear on the original.

EXCEPT AS EXPRESSLY AUTHORIZED ABOVE, CUSTOMER SHALL NOT: COPY, IN WHOLE OR IN PART, SOFTWARE OR DOCUMENTATION; MODIFY THE SOFTWARE; REVERSE COMPILE OR REVERSE ASSEMBLE ALL OR ANY PORTION OF THE SOFTWARE; OR RENT, LEASE, DISTRIBUTE, SELL, OR CREATE DERIVATIVE WORKS OF THE SOFTWARE.

Customer agrees that aspects of the licensed materials, including the specific design and structure of individual programs, constitute trade secrets and/or copyrighted material of Ludlum. Customer agrees not to disclose, provide, or otherwise make available such trade secrets or copyrighted material in any form to any third party without the prior written consent of Ludlum. Customer agrees to implement reasonable security measures to protect such trade secrets and copyrighted material. Title to Software and documentation shall remain solely with Ludlum.

LIMITED WARRANTY. Ludlum warrants that for a period of ninety (90) days from the date of shipment from Ludlum: (i) the media on which the Software is furnished will be

free of defects in materials and workmanship under normal use; and (ii) the Software substantially conforms to its published specifications. Except for the foregoing, the Software is provided AS IS. This limited warranty extends only to Customer as the original licensee. Customer's exclusive remedy and the entire liability of Ludlum and its suppliers under this limited warranty will be, at Ludlum or its service center's option, repair, replacement, or refund of the Software if reported (or, upon request, returned) to the party supplying the Software to Customer. In no event does Ludlum warrant that the Software is error free or that Customer will be able to operate the Software without problems or interruptions.

This warranty does not apply if the software (a) has been altered, except by Ludlum, (b) has not been installed, operated, repaired, or maintained in accordance with instructions supplied by Ludlum, (c) has been subjected to abnormal physical or electrical stress, misuse, negligence, or accident, or (d) is used in ultrahazardous activities.

DISCLAIMER. EXCEPT AS SPECIFIED IN THIS WARRANTY, ALL EXPRESS OR IMPLIED CONDITIONS, REPRESENTATIONS, AND WARRANTIES INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, NONINFRINGEMENT OR ARISING FROM A COURSE OF DEALING, USAGE, OR TRADE PRACTICE, ARE HEREBY EXCLUDED TO THE EXTENT ALLOWED BY APPLICABLE LAW.

IN NO EVENT WILL LUDLUM OR ITS SUPPLIERS BE LIABLE FOR ANY LOST REVENUE, PROFIT, OR DATA, OR FOR SPECIAL, INDIRECT, CONSEQUENTIAL, INCIDENTAL, OR PUNITIVE DAMAGES HOWEVER CAUSED AND REGARDLESS OF THE THEORY OF LIABILITY ARISING OUT OF THE USE OF OR INABILITY TO USE THE SOFTWARE EVEN IF LUDLUM OR ITS SUPPLIERS HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. In no event shall Ludlum's or its suppliers' liability to Customer, whether in contract, tort (including negligence), or otherwise, exceed the price paid by Customer. The foregoing limitations shall apply even if the above-stated warranty fails of its essential purpose. SOME STATES DO NOT ALLOW LIMITATION OR EXCLUSION OF LIABILITY FOR CONSEQUENTIAL OR INCIDENTAL DAMAGES.

The above warranty DOES NOT apply to any beta software, any software made available for testing or demonstration purposes, any temporary software modules or any software for which Ludlum does not receive a license fee. All such software products are provided AS IS without any warranty whatsoever.


This License is effective until terminated. Customer may terminate this License at any time by destroying all copies of Software including any documentation. This License will terminate immediately without notice from Ludlum if Customer fails to comply with any provision of this License. Upon termination, Customer must destroy all copies of Software.

Software, including technical data, is subject to U.S. export control laws, including the U.S. Export Administration Act and its associated regulations, and may be subject to export or import regulations in other countries. Customer agrees to comply strictly with

all such regulations and acknowledges that it has the responsibility to obtain licenses to export, re-export, or import Software.

This License shall be governed by and construed in accordance with the laws of the State of Texas, United States of America, as if performed wholly within the state and without giving effect to the principles of conflict of law. If any portion hereof is found to be void or unenforceable, the remaining provisions of this License shall remain in full force and effect. This License constitutes the entire License between the parties with respect to the use of the Software.

Restricted Rights - Ludlum's software is provided to non-DOD agencies with RESTRICTED RIGHTS and its supporting documentation is provided with LIMITED RIGHTS. Use, duplication, or disclosure by the Government is subject to the restrictions as set forth in subparagraph "C" of the Commercial Computer Software - Restricted Rights clause at FAR 52.227-19. In the event the sale is to a DOD agency, the government's rights in software, supporting documentation, and technical data are governed by the restrictions in the Technical Data Commercial Items clause at DFARS 252.227-7015 and DFARS 227.7202. Manufacturer is Ludlum Measurements, Inc. 501 Oak Street Sweetwater, Texas 79556.

A dark gray square containing the word "Section" in a white sans-serif font at the top, and a large white number "8" in the center.

Technical Theory of Operation

Refer to Amplifier/Power Supply Board Drawing 390 × 400 for the following:

Detector Input/Amplifier

Negative-going detector pulses are coupled from the detector through C022 to Amplifier U021. R023 and CR021 protect the input of U021 from inadvertent shorts. Self-biased amplifier U021 provides gain in proportion to R021 divided by R022. Transistor pins 4, 5, and 6 of U021, provide amplification. Pins 12 and 15 of U021 are coupled as a constant current source to pin 6 of U021. The output self biases to $2V_{be}$ (approximately 1.4 volts) at pin 7 of U021. This provides just enough bias current through pin 6 of U021 to conduct all of the current for the constant current source. Positive pulses from pin 7 of U021 are coupled to the discriminators through R011 and C011.

Alpha/Beta Discriminator

Positive pulses from amplifier U021 are coupled to comparator U012, pin 6, for alpha discrimination and pins 6 and 2 of U011 for beta discrimination. R201, alpha threshold, provides the reference voltage for alpha comparator U012. R102, beta threshold (defined as the lower threshold limit of the beta counting window), provides the reference voltage for beta threshold comparator pins 1, 2, and 3 of U011. R104, beta window (defined as the upper threshold limit of the beta counting window), provides the reference voltage for the beta window comparator pins 5, 6, and 7 of U011.

Alpha/Beta Discriminator Logic Circuit

Alpha pulses from U012 are coupled to univibrator U111. Pulses at pin 6 of U111 are inverted by Q111 for connection to reset (R) pins 3, 13 of U101. Pin 9 of U111 provides the pulses to be counted by the microprocessor (μP). Pulses from pin 9 of U111 are connected to pin 3 of U111 to provide a time delay for the μP clock cycle to complete before the next alpha pulse can be recognized by the μP .

Beta pulses from pin 1 of U011 are coupled to univibrator U101. Pulses are coupled to the μ P from pin 7 of U101 as long as pins 3 and 13 of U011 remain high (+5 V). When an alpha and/or a beta window pulse is present, the reset (pins 3 and 13 of U101) function is enabled, and 7 of U101 remains high. Pin 7 of U101 is connected to pin 13 of U101 to provide a time delay for the μ P clock cycle to complete before the next beta pulse can be recognized by the μ P.

Low Voltage Supply

Battery voltage is coupled to switching regulator U201 and associated components to provide +5 V to power op-amps and logic circuitry. The charge pump (cp) output C202, CR211, CR212, and C201 form a voltage doubler circuit to provide +9 V for U201 amplifier supply. U001 and related components provide +2.5 V reference for HV SET and alpha/beta discriminator controls. R201 (LO BAT) is adjusted so that the meter pointer is aligned with the left vertical mark on the BAT OK line with 2.2 V battery input.

High Voltage Supply

High voltage is developed by blocking oscillator Q421, T321, C412, and rectified by voltage multiplier CR221-CR224, C221-C223, C211, and C114. High voltage increases as current through Q421 increases, with maximum output voltage with Q421 saturated.

High voltage is coupled back through R123 to op-amp pin 2 of U311. Resistor network R211-214 completes the HV division circuit to ground. R214 provides HV limit at 2.0 kV when the HV SET control on the calibration board is at maximum. The regulated HV output is controlled by HV potentiometer located under the CAL cover on the front panel. This control provides the reference for comparator pin 3, U311. During stable operation, the voltage at pin 2 of U311 will equal the voltage at pin 3 of U311. Pin 1 of U311 will cause conduction of Q312 to increase or decrease until the HV finds a level of stability. R115 (HV TEST) calibrates the analog meter to the HV output when the HV test push-button switch is depressed.

Detector Overload

A voltage drop is developed across R121 and sensed by comparator U012 as detector current increases. When the voltage at pin 3 of U012 goes below pin 2, pin 1 goes low, illuminating the OL LED and driving the meter to full scale. R211, overload, provides adjustment for the overload set point.

Meter Drive

Pulses are coupled from the μ P board to the gate of Q302. Q302 inverts the pulses at CR403, and C401 provides integration. Integrated meter drive voltage is coupled from P1-13 via the battery (BAT) and HV test switch to pin 5 of U311. The meter is driven by the emitter of Q111, coupled as a voltage follower in conjunction with pin 6 and 7 of U311. R406, "Meter Cal," is adjusted to calibrate the ratemeter reading corresponding to the incoming count rate. R407 and R408 provide temperature compensation for changes in the meter resistance due to temperature variations.

Refer to Processor Board Drawing 390 × 395 for the following:

Power supply

Battery voltage is coupled to switching regulator U321 and associated components to provide +5 V to power the μ P and display drivers U211, 212. R101, C101, Q101, and Q201 form a delay switch, which allows U321 to stabilize before the load current is connected to the +5 V supply.

Microprocessor (μ P)

U111, Intel N87C51FC, controls all of the data, control inputs, and display information. The clock frequency is crystal controlled by Y111 and related components at 6.144 MHz. The μ P incorporates internal memory (ROM) storing the program information. C211 resets the μ P at power-up to initiate the start of the program routine. During the program loop the μ P looks at all of the input switches for initiation or status changes and responds accordingly.

The μ P uses pulse-width modulation to control the analog ratemeter. The analog output, RATE (P1-3) is divided into 255 increments in a 166 μ s period. At full meter deflection the low pulse period – leading edge to leading edge – will be 166 μ s, 500 cpm = 163 μ s, 400 cpm = 130 μ s, 200 cpm = 65 μ s, 100 cpm = 33 μ s, and 0 = no pulse or +5 V. The pulses are inverted by Q302 on the amp/power supply board and then integrated by R403, C401.

LCD Drive

U101 and U001 make up the liquid crystal display drive circuitry. The display information is sent from the μ P to U101 and U001 via DATA 0-1 lines. Each bit is latched into the drivers when the CLOCK line is brought high, then low by the μ P. When 32 bits have been clocked to the drivers, the

LOAD line is brought high, and then low. The corresponding digits and segments are illuminated, corresponding to the stored-count information from the μP .

Audio

Alpha and/or beta audio pulse frequency is generated by the μP and coupled to Q211. Q211 then inverts the pulses and drives the low side of the audio transformer T321. Front-panel VOL control provides the bias voltage to the top of T321. Secondary winding of T321 is coupled to unimorph speaker via front-panel audio jack.



Troubleshooting

Occasionally, you may encounter problems with your LMI instrument or detector that may be repaired or resolved in the field, saving turn-around time and expense in returning the instrument to us for repair. Toward that end, LMI electronics technicians offer the following tips for troubleshooting the most common problems. Where several steps are given, perform them in order until the problem is corrected. Keep in mind that the most common problems encountered with this particular instrument are: (1) detector cables, (2) sticky meters, (3) battery contacts.

Note that the first troubleshooting tip is for determining whether the problem is with the electronics or with the detector. A Ludlum Model 500 Pulser can be invaluable at this point because of its ability to simultaneously check high voltage, input sensitivity or threshold, and the electronics for proper counting.

We hope these tips will prove to be helpful. As always, please call if you encounter difficulty in resolving a problem or if you have any questions.

Troubleshooting Electronics that Utilize Proportional and Scintillator Type Detectors

| <u>SYMPTOM</u> | <u>POSSIBLE SOLUTION</u> |
|--|---|
| No power (or meter does not reach BAT TEST or BAT OK mark) | <ol style="list-style-type: none"> 1. Check batteries and replace if weak. 2. Check polarity (see marks inside battery lid). Are the batteries installed backwards? |

| <u>SYMPTOM</u> | <u>POSSIBLE SOLUTION</u> |
|--|--|
| No power (or meter does not reach BAT TEST or BAT OK mark) (continued) | <ol style="list-style-type: none">3. Check battery contacts. Clean them with rough sandpaper or use an engraver to clean the tips.4. Check for loose or broken wires, especially between the main board and the calibration board. |
| Nonlinear Readings | <ol style="list-style-type: none">1. Check the high voltage (HV) by using a Ludlum Model 500 Pulser (or equivalent). If a multimeter is used to check the HV, ensure that one with high impedance is used, as a standard multimeter could be damaged in this process.2. Check for noise in the detector cable by disconnecting the detector, placing the instrument on the lowest range setting, and wiggling the cable while observing the meter face for significant changes in readings.3. Check for “sticky” meter movement. Does the reading change when you tap the meter? Does the meter needle “stick” at any spot?4. Check the “meter zero.” Turn the power OFF. The meter should come to rest on “0.” |
| Meter goes full scale or “pegs out” | <ol style="list-style-type: none">1. Replace the detector cable to determine whether or not the cable has failed, causing excessive noise.2. Check the HV and, if possible, the input threshold for proper setting. |

SYMPTOM

Meter goes full scale
or “pegs out”
(continued)

POSSIBLE SOLUTION

3. Open the instrument can and check for loose wires.
4. Ensure that the instrument’s can is properly attached. When attached properly, the speaker will be located on the left side of the instrument. If the can is on backwards, interference between the speaker and the input preamplifier may cause noise.

Section 10

Recycling

Ludlum Measurements, Inc. supports the recycling of the electronics products it produces for the purpose of protecting the environment and to comply with all regional, national, and international agencies that promote economically and environmentally sustainable recycling systems. To this end, Ludlum Measurements, Inc. strives to supply the consumer of its goods with information regarding reuse and recycling of the many different types of materials used in its products. With many different agencies – public and private – involved in this pursuit, it becomes evident that a myriad of methods can be used in the process of recycling. Therefore, Ludlum Measurements, Inc. does not suggest one particular method over another, but simply desires to inform its consumers of the range of recyclable materials present in its products, so that the user will have flexibility in following all local and federal laws.

The following types of recyclable materials are present in Ludlum Measurements, Inc. electronics products and should be recycled separately. The list is not all-inclusive, nor does it suggest that all materials are present in each piece of equipment:

| | | |
|----------------|----------|------------------------------|
| Batteries | Glass | Aluminum and Stainless Steel |
| Circuit Boards | Plastics | Liquid Crystal Display (LCD) |

Ludlum Measurements, Inc. products, which have been placed on the market after August 13, 2005, have been labeled with a symbol recognized internationally as the “crossed-out wheelie bin.” This notifies the consumer that the product is not to be mixed with unsorted municipal waste when discarding; each material must be separated. The symbol will be placed near the AC receptacle, except for portable equipment where it will be placed on the battery lid.

The symbol appears as such:



Section 11

Parts List

| | <u>Reference</u> | <u>Description</u> | <u>Part Number</u> |
|--|------------------|---|--------------------|
| Model 2360 Scaler/Ratemeter | UNIT | Completely Assembled Model 2360 Scaler/Ratemeter | 48-2872 |
| Amplifier/Power Supply Board, Drawing 390 × 400 | BOARD | Completely Assembled Amplifier/ Power Supply Board | 5390-400 |
| CAPACITORS | C001 | 100pF, 100V | 04-5661 |
| | C002 | 47pF, 100V | 04-5660 |
| | C011 | 0.1μF, 50V | 04-5663 |
| | C012-C013 | 0.001μF, 100V | 04-5659 |
| | C016 | 0.001μF, 100V | 04-5659 |
| | C017 | 10μF, 35V | 04-5655 |
| | C021 | 10pF, 100V | 04-5673 |
| | C022 | 100pF, 3KV | 04-5532 |
| | C1 | 1μF, 16V | 04-5701 |
| | C101-C102 | 47pF, 100V | 04-5660 |
| | C111-C113 | 47pF, 100V | 04-5660 |
| | C114 | 0.0047μF, 3KV | 04-5547 |
| | C121-C122 | 0.0047μF, 3KV | 04-5547 |
| | C201-C202 | 10μF, 35V | 04-5655 |
| | C203 | 330pF, 100V | 04-5657 |
| | C211 | 0.0047μF, 3KV | 04-5547 |
| | C212 | 68μF, 10V | 04-5654 |
| | C213 | 1μF, 35V | 04-5656 |
| | C214 | 0.01μF, 50V | 04-5664 |
| | C221-C223 | 0.0047μF, 3KV | 04-5547 |
| | C301 | 68μF, 10V | 04-5654 |
| | C311 | 0.01μF, 100V | 04-5523 |

| | <u>Reference</u> | <u>Description</u> | <u>Part Number</u> |
|------------------------|------------------|---------------------|--------------------|
| | C401 | 0.1 μ F, 50V | 04-5663 |
| | C411 | 0.1 μ F, 50V | 04-5663 |
| | C412 | 1 μ F, 35V | 04-5656 |
| | C421 | 68 μ F, 10V | 04-5654 |
| TRANSISTORS | Q111 | 2N7002L | 05-5840 |
| | Q301 | MMBT4403LT1 | 05-5842 |
| | Q302 | 2N7002L | 05-5840 |
| | Q311-Q312 | MMBT3904T | 05-5841 |
| | Q421 | MJD210 | 05-5843 |
| INTEGRATED CIRCUITS | U001 | LM285MX-2.5 | 06-6291 |
| | U011-U012 | TLC372ID | 06-6290 |
| | U021 | CMXT3906 | 05-5890 |
| | U022-U023 | CMXT3904 | 05-5888 |
| | U101 | CD74HC4538M | 06-6297 |
| | U111 | CD74HC4538M | 06-6297 |
| | U201 | MAX631AESA | 06-6285 |
| | U301 | CD74HC4066M | 06-6323 |
| | U311 | TLC27M7ID | 06-6292 |
| DIODES | CR021 | MMBD7000LT1 | 07-6355 |
| | CR111-CR112 | MMBD914L | 07-6353 |
| | CR211-CR212 | BAT54 | 07-6354 |
| | CR221-CR225 | GI250-2 | 07-6266 |
| | CR411 | MMBD914L | 07-6353 |
| THERMISTOR | R407 | 250K, 250mW, 1% | 07-6366 |
| POTENTIOMETERS | R102 | 10K, BETA THRESH | 09-6921 |
| | R103 | 1M, BETA WIN | 09-6906 |
| | R112 | 1G, FHV-1, 2% | 12-7686 |
| | R115 | 1M, HV READOUT | 09-6909 |
| | R201 | 1M, ALPHA THRESH | 09-6906 |
| | R202 | 200K, LO BAT (LB) | 09-6908 |
| | R211 | 1M, OVERLOAD | 09-6906 |
| | R214 | 1M, HV LIMIT | 09-6906 |
| | R406 | 5K, METER CAL (MTR) | 09-6907 |

| | <u>Reference</u> | <u>Description</u> | <u>Part Number</u> |
|-----------|------------------|--------------------|--------------------|
| RESISTORS | R001 | 22.1K, 250mW, 1% | 12-7843 |
| | R002 | 249K, 250mW, 1% | 12-7862 |
| | R003 | 22.1K, 250mW, 1% | 12-7843 |
| | R004 | 1.5K, 250mW, 1% | 12-7878 |
| | R011 | 100 Ohm, 250mW, 1% | 12-7840 |
| | R012 | 22.1K, 250mW, 1% | 12-7843 |
| | R013 | 33.2K, 250mW, 1% | 12-7842 |
| | R014 | 10.0K, 250mW, 1% | 12-7839 |
| | R015 | 22.1K, 250mW, 1% | 12-7843 |
| | R016 | 10.0K, 250mW, 1% | 12-7839 |
| | R021 | 392K, 250mW, 1% | 12-7841 |
| | R022-R023 | 10.0K, 250mW, 1% | 12-7839 |
| | R024 | 33.2K, 250mW, 1% | 12-7842 |
| | R025 | 22.1K, 250mW, 1% | 12-7843 |
| | R026 | 1.00M, 250mW, 5% | 10-7028 |
| | R101 | 100K, 250mW, 1% | 12-7834 |
| | R104 | 22.1K, 250mW, 1% | 12-7843 |
| | R105 | 100K, 250mW, 1% | 12-7834 |
| | R111 | 100 Ohm, 250mW, 1% | 12-7840 |
| | R112 | 1G, FHV-1, 2% | 12-7686 |
| | R113-R114 | 100K, 250mW, 1% | 12-7834 |
| | R116 | 249K, 250mW, 1% | 12-7862 |
| | R121 | 4.7M, 250mW, 5% | 10-7030 |
| | R122 | 1.00M, 250mW, 5% | 10-7028 |
| | R123 | 1G, FHV-1, 2% | 12-7686 |
| | R212-R213 | 1.00M, 250mW, 1% | 12-7844 |
| | R215 | 1.00M, 250mW, 1% | 12-7844 |
| | R301 | 2.21K, 250mW, 1% | 12-7835 |
| | R302-R303 | 22.1K, 250mW, 1% | 12-7843 |
| | R311 | 10.0K, 250mW, 1% | 12-7839 |
| | R312 | 22.1K, 250mW, 1% | 12-7843 |
| | R313 | 2.21K, 250mW, 1% | 12-7835 |
| | R314 | 10.0K, 250mW, 1% | 12-7839 |
| | R401 | 200 Ohm, 250mW, 1% | 12-7846 |
| | R402 | 221K, 250mW, 1% | 12-7845 |
| | R403 | 7.5K, 250mW, 1% | 12-7847 |
| | R404 | 2.21K, 250mW, 1% | 12-7835 |
| | R405 | 1.00M, 250mW, 1% | 12-7844 |
| | R408 | 301 Ohm, 250mW, 1% | 12-7863 |
| | R411 | 200 Ohm, 250mW, 1% | 12-7846 |
| | R412 | 10.0K, 250mW, 1% | 12-7839 |

| | <u>Reference</u> | <u>Description</u> | <u>Part Number</u> |
|---|------------------|---|--------------------|
| CONNECTORS | P1 | 1-640456-5 MTA 100x15 | 13-8355 |
| | P2 | 640456-3 MTA 100x3 | 13-8081 |
| INDUCTOR | L301 | 220 μ H | 21-9678 |
| TRANSFORMER | T321 | L8050 x50 | 40-0902 |
| Processor Board, Drawing 390 x 395 | BOARD | Completely Assembled Processor Board | 5390-395 |
| CRYSTALS | Y111 | MICRO XTAL-6.144 MHZ | 01-5262 |
| | Y401 | MICRO XTAL-32.768 KHZ | 01-5305 |
| CAPACITORS | C111-C112 | 27PF, 100V | 04-5658 |
| | C211 | 4.7 μ F, 25V | 04-5653 |
| | C221 | 47 μ F, 16V | 04-5666 |
| | C311-C316 | 68 μ F, 10V | 04-5654 |
| | C401 | 10PF, 100V | 04-5673 |
| | C411 | 10 μ F, 35V | 04-5655 |
| | C421 | 4.7 μ F, 25V | 04-5653 |
| | C422 | 10 μ F, 35V | 04-5655 |
| | C423 | 4.7 μ F, 25V | 04-5653 |
| TRANSISTORS | Q101 | TRANS-2N7002L | 05-5840 |
| | Q102 | TRANS-2N7002L | 05-5840 |
| | Q211 | TRANS-2N7002L | 05-5840 |
| INTEGRATED CIRCUITS | U6 | LT1304CS8 | 06-6394 |
| | U101 | 24C65T-I | 06-6401 |
| | U111 | AT89C51RC2 | 06-6893 |
| | U201 | 24C65T-I | 06-6401 |
| | U202 | PCF8574TD | 06-6402 |
| | U401 | PCF8593TD | 06-6403 |
| | U421 | MAX220ESE | 06-6329 |

| | <u>Reference</u> | <u>Description</u> | <u>Part Number</u> |
|-------------|-------------------|----------------------|--------------------|
| DIODES | CR211-CR212 | CXSH-4LF | 07-6358 |
| | CR401-CR402 | MMBD914L | 07-6353 |
| SWITCH | S301 | CONFIGURE SWITCH | 08-6710 |
| RESISTORS | R101-R102 | 200 Ohm, 125mW, 1% | 12-7846 |
| | R103 | 22.1K, 250mW, 1% | 12-7843 |
| | R104 | 100K, 250mW, 1% | 12-7834 |
| | R201 | 22.1K, 250mW, 1% | 12-7843 |
| | R202 | 100K, 250mW, 1% | 12-7834 |
| | R211 | 100K, 250mW, 1% | 12-7834 |
| | R311 | 150K, 250mW, 1% | 12-7833 |
| | R312 | 22.1K, 250mW, 1% | 12-7843 |
| | R313 | 68.1K, 250mW, 1% | 12-7881 |
| | R314 | 1M, 250mW, 1% | 12-7844 |
| | R321 | 0 Ohm, 250mW, 1% | 12-7104 |
| | R401 | 1.00K, 250mW, 1% | 12-7832 |
| | RESISTOR NETWORKS | RN201-RN202 | 220K, EXB2HV224JV |
| CONNECTORS | P3 | 1-640456-7 MTA100x17 | 13-8121 |
| | P4 | 1-640456-1 MTA100x11 | 13-8059 |
| | P12 | 640456-6 MTA100x6 | 13-8095 |
| | P13 | 640456-2 MTA100x2 | 13-8073 \ |
| INDUCTOR | L211 | 22 μ H | 21-9808 |
| BATTERY | B411 | 3V LITHIUM | 22-9794 |
| TRANSFORMER | T321 | XFMR- M 177 AUDIO | 4275-083 |

| | <u>Reference</u> | <u>Description</u> | <u>Part Number</u> |
|---|------------------|---|--------------------|
| LCD Display Board, Drawing 390 × 372 | BOARD | Completely Assembled LCD Display Board | 5390-372 |
| CAPACITOR | C201 | 27PF, 100V | 04-5658 |
| INTEGRATED CIRCUITS | U001 | AY0438-I/L | 06-6358 |
| | U101 | AY0438-I/L | 06-6358 |
| LEDS | DS001 | E118, RED | 07-6308 |
| | DS201 | E118, RED | 07-6308 |
| LCD | DSP101 | LCD-GD-7427WP 6 DIGIT | 07-6624 |
| RESISTORS | R011 | 10K, 125mW, 1% | 12-7839 |
| | R012 | 10K, 125mW, 1% | 12-7839 |
| | R111 | 10K, 125mW, 1% | 12-7839 |
| | R112 | 10K, 125mW, 1% | 12-7839 |
| Calibration Board, Drawing 390 × 176 | BOARD | Completely Assembled Calibration Board | 5390-172 |
| POTENTIOMETER | R3 | 250K, HV SET | 09-6819 |
| CONNECTOR | P7 | 640456-5 MTA100 | 13-8057 |
| BCD Board, Drawing 261 × 107 | BOARD | Completely Assembled BCD Board | 5261-154 |
| SWITCH | S111 | 513384 | 08-6656 |
| CONNECTOR | P5 | 640456-5 MTA100 | 13-8057 |

Chassis Wiring
Diagram, Drawing
390 x 179

| | <u>Reference</u> | <u>Description</u> | <u>Part Number</u> | |
|---------------|------------------|-------------------------|----------------------|---------|
| SWITCHES | S1 | SWITCH-PA-600-210 | 08-6501 | |
| | S3 | SWITCH-MPS-103F | 08-6699 | |
| | S4 | PHONE JACK TINI #42A | 21-9333 | |
| | S5 | SWITCH-7103SYZQE TOGGLE | 08-6720 | |
| | S6 | SWITCH-7205SYZQE TOGGLE | 08-6750 | |
| | POTENTIOMETER | R1 | 10K, VOLUME | 09-6753 |
| CONNECTORS | J1 | CONN-1-640442-5 MTA100 | 13-8383 | |
| | J2 | CONN-640442-3 MTA100 | 13-8135 | |
| | J3 | CONN-1-640442-7 MTA100 | 13-8505 | |
| | J4 | CONN-640442-3 MTA100 | 13-8135 | |
| | J5 | CONN-640442-6 MTA100 | 13-8171 | |
| | J6 | D RECPT-RD9F000V3 9 PIN | 13-8003 | |
| | J7 | JACK-09-9011-1-0419 | 18-9080 | |
| | J8 | CONN-640442-6 MTA100 | 13-8171 | |
| | J9 | CONN-640442-2 MTA100 | 13-8178 | |
| | J12 | CONN-640442-6 MTA100 | 13-8171 | |
| | J13 | CONN-640442-2 MTA100 | 13-8178 | |
| | J14 | Series "C"-UG706/U | 13-7751 | |
| | BATTERIES | B1-B2 | "D" DURACELL BATTERY | 21-9313 |
| | AUDIO | DS1 | S100RL-M, UNIMORPH | 21-9676 |
| MISCELLANEOUS | M1 | METER ASSEMBLY | 4390-160 | |
| | * | METER-PORT BZL W/GLS | 4363-352 | |
| | * | M2360 METERFACE | 7390-157 | |
| | * | METER-MVT #919492 1 MA | 15-8030 | |
| | * | O RING-BEZEL | 16-8334 | |
| | * | SPACER-#2 X .187 NYL | 18-9143 | |
| | * | BATTERY CONTACT SET | 40-1707 | |
| | * | MAIN HARNESS | 8390-162 | |
| | * | LCD HARNESS | 8390-163 | |
| | * | CAN-RS232 HARNESS | 8390-166 | |

A dark gray square graphic with the word "Section" in a bold, sans-serif font at the top. Below it, the number "12" is written in a very large, bold, white sans-serif font.

Drawings and Diagrams

AMPLIFIER/POWER SUPPLY BOARD, Drawing 390 × 400

AMPLIFIER/POWER SUPPLY BOARD COMPONENT LAYOUT,
Drawing 390 × 403 (2 sheets)

PROCESSOR BOARD, Drawing 390 × 395

PROCESSOR BOARD COMPONENT LAYOUT,
Drawing 390 × 398 (2 sheets)

LCD DISPLAY BOARD, Drawing 390 × 372

LCD DISPLAY BOARD COMPONENT LAYOUT, Drawing 390 ×
373A (2 sheets)


CALIBRATION BOARD, Drawing 390 × 176

CALIBRATION BOARD COMPONENT LAYOUT,
Drawing 390 × 177 (2 sheets)

BCD BOARD, Drawing 261 × 107

BCD BOARD COMPONENT LAYOUT, Drawing 261 × 105A
(2 sheets)


WIRING DIAGRAM, Drawing 390 × 179



| | | | | | |
|--|--|---|--|--|------------|
| | | LUDLUM MEASUREMENTS, INC. | | PO Box 816 801 Oak Street Sweetwater, Texas 79556 U.S.A. 1-800-822-0828 | |
| | | Drawn: RDS 11/20/2019 Design: DAW 11/20/2019 Approved: <i>[Signature]</i> 11/20/2019 Print Date: 11/21/2019 9:27:27 AM | Title: AMPLIFIER/PS BOARD Model: 2360 Board: S390-400 Sheet: 1 of 1 Rev: 1 | Series: 390 | Sheet: 400 |

DO NOT PLACE C500 AND R500

5390-400
REV: 1



LUDLUM MEASUREMENTS, INC.
 PO Box 810
 501 Oak Street
 Sweetwater, TX 79556
 U.S.A. 1-800-622-0828


Title: AMPLIFIER/PS BOARD

| | | |
|-----------------------------------|------------------|----------------------|
| Drawn: PAB | 11/20/2019 | Model: 2360 |
| Design: RDS | 11/20/2019 | Board#: 5390-400 |
| Approve: <i>RYS</i> | <i>21 Nov 19</i> | Rev: 1 |
| PCBA Drawing | | SCALE: 1.08 |
| Print Date: 11/21/2019 9:27:34 AM | Top Overlay | Series 390 Sheet 403 |

\\freedom\pc\Projects\LMI\5390-400\Rev1\390400R1_Assy_PcbDoc

5390-400

REV: 1



+




LUDLUM
MEASUREMENTS, INC.

PO Box 810
501 Oak Street
Sweetwater, TX 79556
U.S.A. 1-800-622-0828

| | | | |
|-----------------------------------|------------|------------------|--------------|
| Title: AMPLIFIER/PS BOARD | | | |
| Drawn: PAB | 11/20/2019 | Model: 2360 | |
| Design: RDS | 11/20/2019 | Board#: 5390-400 | |
| Approve: <i>RJS 21NOV19</i> | | Rev: 1 | |
| PCBA Drawing | | SCALE: 1.08 | Series Sheet |
| Print Date: 11/21/2019 9:27:38 AM | | Bottom Overlay | 390 403 |

\\freedom\pcb\Projects\LMMIM 2360\5390-400\Rev1\09040R1_Assy.PcbDoc




LUDLUM MEASUREMENTS, INC.

PO Box 810
801 Oak Street
Sweetwater, Texas 79556
U.S.A. 1-800-822-0823

| | | |
|------------------|------------|------------------------|
| Drawn: PAB | 11/22/2019 | Title: PROCESSOR BOARD |
| Design: RSS | 11/22/2019 | Model: 2360 |
| Board#: 5390-395 | | |
| Sheet: 1 of 1 | | |
| Series | | Sheet |
| 390 | | 395 |

Print Date: 11/25/2019 1:24:44 PM Rev: 1



5390-395
REV: 1

LUDLUM MEASUREMENTS, INC. PO Box 810
501 Oak Street
Sweetwater, TX 79556
U.S.A. 1-800-622-0828


Title: PROCESSOR BOARD

| | | |
|-----------------------------|--------------------|------------------|
| Drawn: PAB | 11/22/2019 | Model: 2360 |
| Design: RSS | 11/22/2019 | Board#: 5390-395 |
| Approve: <i>[Signature]</i> | <i>[Signature]</i> | Rev: 1 |

PCBA Drawing

| | | | |
|-----------------------------------|-------------|------------|-----------|
| Print Date: 11/25/2019 1:26:31 PM | SCALE: 1.08 | Series 390 | Sheet 398 |
| Top Overlay | | | |

\\freedom\pc\Projects\LMM\2360\5390-395\Rev1\390395R1_Assy.PcbDoc



5390-395
REV: 1

| | | | |
|--|-----------------|---|--------------|
|  LUDLUM MEASUREMENTS, INC. | | PO Box 810 501 Oak Street Sweetwater, TX 79556 U.S.A. 1-800-622-0828 | |
| Title: PROCESSOR BOARD | | | |
| Drawn: PAB | 11/22/2019 | Model: 2360 | |
| Design: RSS | 11/22/2019 | Board#: 5390-395 | |
| Approve: <i>PAB</i> | <i>11/22/19</i> | Rev: 1 | |
| PCBA Drawing | | SCALE: 1.08 | Series Sheet |
| Print Date: 11/25/2019 1:26:31 PM | | Bottom Overlay | 390 398 |
| <small>\\freedom\pcb\Projects\LM\IM 2360\5390-395\Rev1\390395R1_Assembly.PcbDoc</small> | | | |

A

B

C

D

E

1








2

3

4

5


DSP101
S-TEK GD-7427WP




LUDLUM MEASUREMENTS, INC.
 PO Box 810
 501 Oak Street
 Sweetwater, Texas 79556
 U.S.A. 1-800-622-0828

| | | |
|------------------------------|-------------|-------------------------------------|
| Drawn: CKB | 11/08/2013 | Title: LCD Display -STEK LCD (2013) |
| Design: RSS | 11/08/2013 | Model: 2360 |
| Approved: <i>[Signature]</i> | 11/08/2013 | Board#: 5390-372 |
| Print Date: 7/1/2014 | 11:30:52 AM | Rev: 2 |
| Series | Sheet | |
| 390 | 372 | |

W:\Projects\LMM 2360\5390-372\Rev2\390372R2P1 SchDoc




| | | | |
|--|-------------|---|--------------------------|
|  | | PO Box 810 501 Oak Street Sweetwater, TX 79556 U.S.A. 1-800-622-0828 | |
| Title: LCD Display -STEK LCD (2013) | | | |
| Drawn: CKB | 11/08/2013 | Model: 2360 | |
| Design: RSS | 11/08/2013 | Board#: 5390-372 | |
| Approve: <i>RDS 1 Jul 14</i> | | Rev: 2 | |
| PCBA Drawing | | SCALE: 1.05 | Series |
| Print Date: 7/1/2014 | 11:31:01 AM | Top Overlay | Sheet 390 373A |
| <small>W:\Projects\LMIM 2360\5390-372\Rev2\390372R2_Assy.PcbDoc</small> | | | |




| | | | |
|---|-----------------|---|---------------------|
|  LUDLUM MEASUREMENTS, INC. | | PO Box 810 501 Oak Street Sweetwater, TX 79556 U.S.A. 1-800-622-0828 | |
| | | Title: LCD Display -STEK LCD (2013) | |
| Drawn: CKB | 11/08/2013 | Model: 2360 | |
| Design: RSS | 11/08/2013 | Board#: 5390-372 | |
| Approve: <i>RDS</i> | <i>1 JUL 14</i> | Rev: 2 | |
| PCBA Drawing | | SCALE: 1.05 | Series Sheet |
| Print Date: 7/1/2014 | 11:31:03 AM | Bottom Overlay | 390 373A |
| <small>W:\Projects\LMIM 2360\5390-372\Rev2\390372R2_Assy.PcbDoc</small> | | | |


| REVISIONS | | | | | | |
|-----------|-----------|------|------|-------------|------|----------|
| EFF | AUTHORITY | ZONE | TLTR | DESCRIPTION | DATE | APPROVED |



| | | | | | |
|-------------------|----------------|--------------------------|-------|--------|--------------|
| UPDATED - | 1 / | LUDLUM MEASUREMENTS INC. | | | |
| DR CKB | 12/14/96 | TITLE: CALIBRATION BOARD | | | |
| CHK <i>RSS</i> | <i>3/14/97</i> | BOARD# 5390-172 | | | |
| DSGN RSS | 12/14/96 | SIZE | MODEL | SERIES | SHEET |
| APPD <i>JGW</i> | <i>3/17/97</i> | D | 2360 | 390 | 176 |
| NEXT HIGHER ASSY. | | | | | |
| 13:49:27 | 16-Dec-96 | SB390172 | | | SHEET 1 OF 1 |




| | | |
|--|---------------|---|
|  LUDLUM MEASUREMENTS INC. SWEETWATER, TX. | | |
| DR | CKB 12/14/96 | TITLE: CALIBRATION BOARD |
| | | BOARD# 5390-172 BS390104 |
| DSCN | RSS 12/14/96 | MODEL 2224-1 SERIES 390 SHEET 177 |
| APP | <i>SKG/IA</i> | COMP ARTWORK <input type="checkbox"/> SLDR ARTWORK <input type="checkbox"/> |
| 09:38:12 | 20-Nov-96 | COMP OUTLINE <input type="checkbox"/> SLDR OUTLINE <input checked="" type="checkbox"/> |
| | | COMP PASTE <input type="checkbox"/> COMP MASK <input type="checkbox"/> SLDR PASTE <input type="checkbox"/> SLDR MASK <input type="checkbox"/> |




| | | | | |
|--|------------------|----------|--------------------------|--------------|
|  LUDLUM MEASUREMENTS INC. SWEETWATER, TX. | | | | |
| DR | CKB | 12/14/96 | TITLE: CALIBRATION BOARD | |
| | | | BOARD: 5390-172 | |
| DSCN | RSS | 12/14/96 | MODEL: 2360 | |
| APP | <i>BS/3/1/12</i> | | FILENAME: BS390172 | |
| COMPONENT | SOLDER | | 09:48:56 | 31-JUL-12 |
| | | | REVISION | SERIES SHEET |
| OUTLINE | OUTLINE | | 1.0 | 390 177 |

| POS FUNCTION | |
|--------------|---------|
| 1 | 0.1 MIN |
| 2 | 0.5 MIN |
| 3 | 1 MIN |
| 4 | 2 MIN |
| 5 | 5 MIN |
| 6 | 10 MIN |
| 7 | 60 MIN |
| 8 | PC |
| 10 | |
| 11 | |
| 12 | |



| | | | | |
|---|--|--|----------------------|--|
| | | PO Box 810 501 Oak Street Sweetwater, Texas 79556 U.S.A. 1-800-622-0828 | | |
| | | Drawn: CKB Design: PW | 03/27/01 10/30/96 | Title: BCD Board Model: 2360 Board#: 5261-154 |
| Approve: <i>ROS 7 Sep 12</i> | | Sheet: 1 of 1 | Series | Sheet |
| Print Date: 9/7/2012 4:39:16 PM | | Rev: 1 | 261 | 107 |
| <small>W:\Projects\LMM 2360\5261-154\Rev2.0\261154R2P1.SchDoc</small> | | | | |




LUDLUM MEASUREMENTS, INC.

PO Box 810
501 Oak Street
Sweetwater, TX 79556
U.S.A. 1-800-622-0828

Title: BCD Board

| | | | | |
|--|-----------|-------------------------|---------------|--------------|
| Drawn: CKB | 03/27/01 | Model: 2360 | | |
| Design: PW | 10/30/96 | Board#: 5261-154 | | |
| Approve: RDS | 18 Nov 15 | Rev: 2 | | |
| Print Date: 1/14/2014 2:08:01 PM | | SCALE: 1.00 | Series | Sheet |
| | | Top Overlay | 261 | 105 A |



PX
640456-5



LUDLUM
MEASUREMENTS, INC.
















PO Box 810
501 Oak Street
Sweetwater, TX 79556
U.S.A. 1-800-622-0828

Title: BCD Board

| | | | | |
|--------------------|------------|--------------------|---------------|--------------|
| Drawn: CKB | 03/27/01 | Model: 2360 | | |
| Design: PW | 10/30/96 | Board#: 5261-154 | | |
| Approve: RDS | 18 Nov 15 | Rev: 2 | | |
| Print Date: | | SCALE: 1.00 | Series | Sheet |
| 1/14/2014 | 2:08:04 PM | Bottom Overlay | 261 | 105 A |

W:\Projects\LM\LM 2360\5261-154\Rev2.0\261154R2_Assy.PcbDoc

| REVISONS | | | | | | |
|----------|-----------|------|-----|-------------|------|----------|
| EFF | AUTHORITY | ZONE | LTR | DESCRIPTION | DATE | APPROVED |



| | | | | | |
|-----------------|----------------|--------------------------|--------------|------------|-----------|
| UPDATED | - | LUDLUM MEASUREMENTS INC. | | | |
| DR CKB | 12/20/96 | TITLE: WIRING DIAGRAM | | | |
| CHK <i>RSS</i> | <i>3-19-97</i> | BOARD# 390-173 | | | |
| DSGN RSS | 10/03/96 | SIZE C | MODEL 2360 | SERIES 390 | SHEET 179 |
| APPD <i>Jow</i> | <i>3-19-97</i> | NEXT HIGHER ASSY. | | | |
| 17:27:27 | 10-Feb-97 | W390173 | SHEET 1 OF 1 | | |